

# Technology

## Digestion of Graphite

**Summary** Graphite is digested in a Berghof stainless steel high pressure digestion system.

### Method

#### Equipment

	Type	Manufacturer
	DAB-2	Berghof Products + Instruments GmbH
	Vessel volume 50 mL	Berghof Products + Instruments GmbH

#### Reagents

	Acid	Volume
	HNO <sub>3</sub> (65%)	5.0 mL
	H <sub>2</sub> SO <sub>4</sub> (95%)	5.0 mL

#### Procedure\*

Weigh 500 mg of the sample into the digestion vessel. Add 5.0 mL of HNO<sub>3</sub> and 5.0 mL of H<sub>2</sub>SO<sub>4</sub>. Allow the mixture to react for a period of at least 5 minutes. Close the vessel and heat with the following program.

Temperature Program	Step	T [°C]	Time [h]**
	1	250	48-72

**Note:** To avoid foaming and splashing wait until the vessels have cooled to room temperature. Carefully open the digestion vessel in a fume hood wearing hand, eye and body protection since a large amount of gas will be produced during the digestion process.

#### Results

Clear solution.

**Note:** This application serves only as a guide line and may need to be optimized for your sample.

\*For other high pressure digestion vessels the maximum sample mass and volume of the acid mixture required can easily be calculated from the proportion of the inner vessel volumes.

\*\*This time includes the heating-up period

