

# Technology, Pharma

## Screening of metathesis catalysts

### Summary

Berghof reactor systems are often used for screening of metathesis catalysts. By the use of the reactors it can be tested which effects different catalyst structures have on the catalysis process or which effect smaller changes of process parameters have.

### Examples for Equipment

Equipment	Reactor	Heating instrument	Berghof controller	Stirring via hot plate available
	BR-25	BLH/BAH	BTC-3000	✓
	BR-40	BLH/BAH	BTC-3000	✓
	BR-100	BLH/BAH	BTC-3000	✓
	BR-200	BLH/BAH	BTC-3000	✓