



# SAFETY DATA SHEET

according to 29 CFR 1910.1200 and GHS

Page 1/7

## AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product name	AccuPoint® Advanced Alkaline Phosphatase
Product code	9950

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Product Use	[SU3] Industrial uses: Uses of substances as such or in preparations at industrial sites; [PC21] Laboratory chemicals;
Description	Detection test.

#### 1.3. Details of the supplier of the safety data sheet

Company	Neogen Corporation
Address	620 Leshler Place Lansing MI 48912 USA
Web	www.neogen.com
Telephone	517-372-9200/800-234-5333
Email	regulatorysupport@neogen.com

#### 1.4. Emergency telephone number

Emergency telephone number	1-800-234-5333
----------------------------	----------------

#### Further information

Hotline(s):	Medical: United States (1-800-498-5743), International (1-651-523-0318). Spill/CHEMTREC: United States (1-800-424-9300), International (1-703-527-3887).
-------------	--

### SECTION 2: Hazards identification

#### 2.2. Label elements

Hazard Statements	No Significant Hazard
-------------------	-----------------------

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

#### EC 1272/2008

Chemical Name	Index No.	CAS No.	EC No.	REACH Registration Number	Conc. (%w/w)	Classification
Sucrose		57-50-1	200-334-9		70 - 80%	
Magnesium nitrate		10377-60-3	233-826-7		0 - 0.5%	Ox. Sol. 2: H272; Skin Irrit. 2: H315; Eye Irrit. 2: H319; STOT SE 3: H335;

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

## Further information

The exact percentage (concentration) of composition has been withheld as a trade secret.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

Inhalation	Move the exposed person to fresh air.
Eye contact	Rinse immediately with plenty of water for 15 minutes holding the eyelids open.
Skin contact	Wash off immediately with plenty of soap and water. Remove contaminated clothing.
Ingestion	DO NOT INDUCE VOMITING.

### 4.2. Most important symptoms and effects, both acute and delayed

Inhalation	Inhalation may cause coughing, tightness of the chest and irritation of the respiratory system.
Eye contact	May cause irritation to eyes.
Skin contact	May cause irritation to skin.
Ingestion	Ingestion may cause nausea and vomiting.

### 4.3. Indication of any immediate medical attention and special treatment needed

Inhalation	If you feel unwell, seek medical advice (show the label where possible).
Eye contact	Seek medical attention.
Skin contact	Seek medical attention if irritation or symptoms persist.
Ingestion	If you feel unwell, seek medical advice (show the label where possible).

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Use extinguishing media appropriate to the surrounding fire conditions.

### 5.2. Special hazards arising from the substance or mixture

Burning produces irritating, toxic and obnoxious fumes.

### 5.3. Advice for firefighters

Wear self contained breathing apparatus and protective clothing.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation of the working area. Avoid contact with skin and eyes.

### 6.2. Environmental precautions

Do not allow product to enter drains. Prevent further spillage if safe.

### 6.3. Methods and material for containment and cleaning up

Transfer to suitable, labelled containers for disposal. Clean spillage area thoroughly with plenty of water.

### 6.4. Reference to other sections

See section 2, 8, & 13 for further information.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Ensure adequate ventilation of the working area. Avoid contact with eyes and skin. Adopt best Manual Handling considerations when handling, carrying and dispensing.

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

## 7.2. Conditions for safe storage, including any incompatibilities

Keep containers tightly closed. Keep in a cool, dry, well ventilated area. Store in correctly labelled containers.

## 7.3. Specific end use(s)

See section 1.2 for further information.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

Occupational exposure controls.

### Exposure Limit Values

Sucrose	ACGIH TLV TWA: 10 mg/m <sup>3</sup>	ACGIH TLV STEL: -
	Notations: A4	OSHA PEL 8hr TWA ppm: -
	OSHA PEL 8hr TWA mg m3: 15	Skin Designation: -
	NIOSH REL 10hr TWA ppm:	NIOSH REL 10hr TWA mg m3:
	NIOSH REL 15 min ST ppm:	NIOSH REL 15 min ST mg m3:

### 8.2. Exposure controls



Appropriate engineering controls

Ensure adequate ventilation of the working area.

Individual protection measures

Wear protective clothing.

Eye / face protection

Approved safety goggles.

Skin protection - Hand protection

Wear suitable gloves.

Skin protection - Other

Wear suitable protective clothing.

Respiratory protection

Wear suitable respiratory equipment when necessary.

Occupational exposure controls

Keep away from food, drink and animal feed and/ or feed products.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

## 9.1. Information on basic physical and chemical properties

<b>State</b>	Powder
<b>Color</b>	White
<b>Odor</b>	Odorless
<b>Odour threshold</b>	No data available
<b>pH</b>	No data available
<b>Melting point</b>	No data available
<b>Freezing Point</b>	No data available
<b>Boiling point</b>	No data available
<b>Flash point</b>	No data available
<b>Evaporation rate</b>	No data available
<b>Flammability limits</b>	No data available
<b>Vapor pressure</b>	No data available
<b>Vapor density</b>	No data available
<b>Relative density</b>	No data available
<b>Fat Solubility</b>	No data available
<b>Partition coefficient</b>	No data available
<b>Autoignition temperature</b>	No data available
<b>Viscosity</b>	No data available
<b>Explosive properties</b>	No data available
<b>Oxidising properties</b>	No data available
<b>Solubility</b>	Soluble in water

## 9.2. Other information

<b>Conductivity</b>	No data available
<b>Surface tension</b>	No data available
<b>Gas group</b>	No data available
<b>Benzene Content</b>	No data available
<b>Lead content</b>	No data available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

	Stable under normal conditions.
--	---------------------------------

### 10.2. Chemical stability

	Stable under normal conditions.
--	---------------------------------

### 10.3. Possibility of hazardous reactions

	No data is available on this product.
--	---------------------------------------

### 10.4. Conditions to avoid

	No data is available on this product.
--	---------------------------------------

### 10.5. Incompatible materials

	No data available.
--	--------------------

### 10.6. Hazardous decomposition products

	No data is available on this product.
--	---------------------------------------

## SECTION 11: Toxicological information

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

## 11.1. Information on toxicological effects

Skin corrosion/irritation	May cause irritation to skin.
Serious eye damage/irritation	May cause irritation to eyes.
Carcinogenicity	Not found to be a potential carcinogen in the International Agency for Research on Cancer (IARC) Monographs, Volumes 1-112. Not listed in the National Toxicology Program (NTP) 13th Report on Carcinogens. Not listed in the OSHA standard 1910.1003 Carcinogens.

### 11.1.4. Toxicological information

	No data available
--	-------------------

## SECTION 12: Ecological information

### 12.1. Toxicity

	No data available
--	-------------------

### 12.2. Persistence and degradability

	No data is available on this product.
--	---------------------------------------

### 12.3. Bioaccumulative potential

	No data is available on this product.
--	---------------------------------------

### Partition coefficient

	AccuPoint® Advanced Alkaline Phosphatase No data available	Sucrose No data available
--	--	---------------------------

### 12.4. Mobility in soil

	No data is available on this product.
--	---------------------------------------

### 12.5. Results of PBT and vPvB assessment

	No data available.
--	--------------------

### 12.6. Other adverse effects

	No data available.
--	--------------------

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

	Dispose of in compliance with all local and national regulations.
--	---

### Disposal methods

	Dispose of in compliance with all local and national regulations.
--	---

### Disposal of packaging

	Dispose of in compliance with all local and national regulations.
--	---

## SECTION 14: Transport information

### 14.1. UN number

	The product is not classified as dangerous for carriage.
--	--

### 14.2. UN proper shipping name

	The product is not classified as dangerous for carriage.
--	--

### 14.3. Transport hazard class(es)

	The product is not classified as dangerous for carriage.
--	--

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

## 14.4. Packing group

The product is not classified as dangerous for carriage.

## 14.5. Environmental hazards

The product is not classified as dangerous for carriage.

## 14.6. Special precautions for user

The product is not classified as dangerous for carriage.

## 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

The product is not classified as dangerous for carriage.

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

<b>Regulations</b>	<p>All chemicals in this product are TSCA listed or excluded from listing, on the United States Environmental Protection Agency Toxic Substances Control Act (TSCA) inventory. TSCA - United States Toxic Substances Control Act Section 8(b) Inventory.</p> <p>Hazard Communication Standard: 29 Code of Federal Regulations 1910.1200: Occupational Safety and Health Standards: Z: Toxic and Hazardous Substances: Hazard Communication. US State Regulations:</p> <p>New Jersey Right-to-Know Hazardous Substance List: Water: Y Magnesium nitrate: Y Mixture of 5-Chloro-2-methyl-4-isothiazolin and 2-Methyl-2H-isothiazol-3-one ( 3:1): Y</p> <p>Pennsylvania Right-to-Know Hazardous Substance List: Water: Y Magnesium nitrate: Y Copper dinitrate: Y, Magnesium nitrate: Y.</p> <p>Massachusetts Right-to-Know Hazardous Substance List: Magnesium nitrate: Y.</p>
--------------------	--

### Inventories

Ingredient	TSCA	MA	SARA313	CALPROP65	CWA307	CAA112	CWA311	SARA302	SARA304
Sucrose (.alpha.-D-Glucopyranoside, .beta.-D-fructofuranosyl)	Y	N	N	N	N	N	N	N	N

## SECTION 16: Other information

### Other information

<b>Text of Hazard Statements in Section 3</b>	<p>Ox. Sol. 2: H272 - May intensify fire; oxidiser.</p> <p>Skin Irrit. 2: H315 - Causes skin irritation.</p> <p>Eye Irrit. 2: H319 - Causes serious eye irritation.</p> <p>STOT SE 3: H335 - May cause respiratory irritation.</p>
---	--

### Further information

	<p><b>DISCLAIMER:</b> The information and recommendations set forth herein ("Information") are presented in good faith and believed to be correct as of the date issued. No representation is made regarding the completeness or accuracy of the Information. Further, because of the many factors that affect the use of this product, the Information is supplied upon the condition that the person(s) receiving it will make their own determination regarding its suitability for their own unique purpose(s), prior to use. Except as expressly stated herein, <b>NO WARRANTIES, GUARANTEES, OR REPRESENTATIONS OF ANY KIND, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, COURSE OF PERFORMANCE, USAGE OF TRADE, OR RESULTS TO BE OBTAINED BY USE OF THIS PRODUCT</b> are made with respect to this product or use of this product. The product covered is furnished "as is" and only subject to the warranties herein provided, no liability is assumed resulting from the use of this product.</p>
--	---

# AccuPoint® Advanced Alkaline Phosphatase

Revision 0  
Revision date 2018-05-18

---