

Option Jigs

別途治具が必要な場合は、お気軽にご相談ください。

D	T
10	12
15	12
20	12
30	12
40	12
50	20
60	30

Flat indenter
Made by POM
CJ-FID-□

θ	D
45	60
45	60
60	30

Cone indenter
Made by POM
CJ-CID-□

θ
45
60

Wedge indenter
Made by POM
CJ-VID-30

Joint Shaft
ジョイントシャフト
GJ-M6-L120

Semi cylinder indenter
POM製
CJ-H15-70

D	T
3	-
5	0.8
8	0.8
10	1
12.7	3
15	3
20	3
30	5
50	5

Flat indenter
Made by SUS
C□-FD□

D
5
7
10
20

Flat indenter
Made by SUS
CT-B□

A/LKB-F
Knife Blade
CT-KB

W
45
80
100

Knife holder
Made by Nylon
CJ-K□

Wire Cutter
Made by SUS
CT-WC

W
5
10
20
30
40

Small Chuck
GS-□

Universal joint
Mde by SUS
GU-M6N6

Noodle stretching test jig
CJ-N20

JIS K6503 Jelly strength test jigs

Flat type 12.7-L35
Made by POM
CJ-F127-35

Jelly cup
CS-G60